EXPANDABLE CORPECTOMY CAGE

GENİŞLEYEBİLİR KORPEKTOMİ KAFES

Introduction

Giriş

EN The Expandable Corpectomy Cages may be used after complete or partial destruetions of the vertebral body. Such destruetions may be caused by; •Tumours •Fractures •infeetions

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EN Expandable Corpectomy Cage is applied if the anterior thoracic or lumbar spine needs to be reconstructed.

Features of the Expandable Corpectomy Cage;

- One or more vertebral bodies can be substituted
- Expandable Corpectomy Cage may be distracted continuously in situ •Expandable Corpectomy Cage corrects malalignements of the spine by distraction •Expandable Corpectomy Cage is anchored easily and safely between the adjacent end plates of the vertebral bodies.
- A small number of instruments makes handling easy.
- The Corpectomy Cage has a hollow design. If necessary the implant may be filled with bone graft through large orifices in the cylindrical portions. This may be done either before or after insertion.

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EN The Expandable Corpectomy Cage consists of three components: a center piece and two end pieces. The implant is available in three external diameters, 20 mm and 24 mm. The center piece's one side is equipped with a right and the otherside is equipped with a left thread. Accordingly, the end pieces are equipped with either a left or a right thread which are marked with "R" and "L".

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Unilateral Expandable Corpectomy Cage

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EN The Unilateral Expandable Corpectomy Cage has already been used for some years as a vertebral body replacement for the anterior cervical and thoracic vertebral column. Unilateral Expandable Corpectomy Cage can also be used for intersomatic fusion.

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EN Unilateral Expandable Corpectomy Cage is used for the reconstruction of substance defects, support of the anterior column as well as restoration of the vertebral column profile. An additional anterior or posterior stabilization by means of a rod/screw system or a plate system is necessary.

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EN By turning the distraction ring the implant is spread open and the distraction height becomes fixed. The distraction height of the Unilateral Corpectomy may be adjusted continuously as per the individual case.

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Material Malzeme

EN Certified Ti6Al4Veli alloy according to ASTM F 136 - ISO 5832-3. It is selected for its, modulus of elasticity closest to that living bone, high fatigue strength and perfect biocompatibility.

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Please Note : This document is intended as a guide for the surgeon only. There are multiple techniques for the insertion of Expandable Corpectomy Cage and, as with any surgical procedure, a surgeon should be thoroughly trained and beware that this procedure is appropriate for the patient before proceeding.

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Enposure Açma İşlemi

EN Determine the correct approach. The involved segment(s) of the vertebral column are exposed. Complete or partial resection of the vertebral body (bodies) including the adjacent discs is performed.

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EN The gap to be bridged is measured using the measuring caliper (Ref: 1068306001).

Corresponding size with measuring calliper on the measuring template (Ref:1083070001) shows the appropriate size of the implant.

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
A Attention / Dikkat

EN The first section on the measuring template corresponds to the Expandable Corpectomy Cage that has not been distracted. The second(striped) section illustrates the possible distraction height.

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EN Unilateral Expandable Corpectomy Cage's measuring table is below the Expandable Corpectomy Cage's table on the measuring template. Chose the appropriate Unilateral Expandable Corpectomy Cage to the measuring caliper.

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
Distraksiyon

EN The expansion instrument (Ref:10683010005) is inserted into the rotation holes of the center piece. Distraction is achieved by turning the center piece as indicated by the arrows.

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
EN As an alternative to the Expansion Instrument the Expansion Instrument/5 may be used. The expansion instrument/5 is attached to the rotation hole of the center piece of the cage. To achieve desired distraction, expansion instrument/5 turned in the direction of the arrow.

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As soon as the appropriate distraction height has been achieved the fixation instrument (Ref: 10683050001) is removed.

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
EN After the correct position of the implant has been checked the end pieces are locked to the center piece using a set screw (Ref:10133202010). Again the screw driver for setscrew (Ref:10683030001) is used.

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

**Bone Grafting
KemikGrafti**

EN The void in the Expandable Corpectomy Cage resulting from the distraction can be topped up through the large drill holes in the end pieces. Additional bone grafting around the implant is possible. The Expandable Corpectomy Cage may also be filled with bone cement depending on the situation.

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
Instruments / El Aletleri

10683010005		Expansion Instrument
10683020002		Expansion Instrument/2
10683020003		Expansion Instrument/3
10683020005		Expansion Instrument/5
10683030001 10683040001		Screw Driver for Setscrew
10683040002		Impactor/Straight
10683050001		Impactor/Angled
10683060001		Fixation Instrument
		Measuring Caliper
10683070001		Measuring Table

**SOLİD TIBBİ MALZEMELER TURİZM İNŞ.
İTH. SAN ve TİC. LTD. ŞTİ**

Adres: Yenidoğan mah. 1525 Sok. no:9 k:1 d:1 54100
Adapazarı /SAKARYA / TÜRKİYE

Tel: +90 264 274 69 63

Tel2: +90 264 278 72 19