

SOLID

a
major
resto-
ration

for the
master-
piece
human





Cover Image

The Caesar Augustus Statue

Gulf of Naples, Island of Ischia, Italy

Kapak görseli

Sezar Augustus Heykeli

Napoli Körfezi, Ischia Adası, İtalya

Изображение на обложке

Статуя императора Августа

Неаполитанский залив, остров Искья, Италия

Bilgi: Bu katalogda göreceğiniz sanat eserleri tasarımın etkisini artırmak amacıyla manipüle edilmiştir. Görsellerin, eserlerin günümüzdeki durumları ile ilgisi yoktur. / **Information:** The artworks available in this catalog have been manipulated to enhance the design effect. The visuals or artworks do not reflect their present status. / **Инфо:** Изображения шедевров искусства, которые вы увидите в этом каталоге, были обработаны с целью повышения качества дизайна. Изображения не связаны с текущим состоянием памятников истории.



34. Uluslararası Liderlik ve Prestij Şirket Ödülleri, "Yeni Milenyum Ödülü", İsviçre

34-я Международная церемония награждения лидирующих и престижных компаний,
«Награда нового тысячелетия», Швейцария



The Great Sphinx of Giza
2723-2563 BC, Cairo, Egypt

Büyük Gize Sfenksi
MÖ 2723-2563, Kahire, Mısır

Большой сфинкс в Гизе
2723-2563 до н.э., Каир, Египет

MP Calcaneal Plate



Product Description

It is indicated for fractures and osteotomies.
Anatomically shaped,
Low Profile Locking plate,
three size variants are key
features of this plate

Screw used:
3.5mm LHS Screw,
3.5mm Flat Head Cortex.

MP Calcaneal Plate

Codes			Range
Small	Medium	Large	
936.L.001	936.L.002	936.L.003	Left
936.R.001	936.R.002	936.R.003	Right

MP Stick Plate

Product Description

Codes

Range

Holes

939.225

22.5mm

02

939.275

27.5mm

02

939.325

32.5mm

02

939.370

37.0mm

03

939.500

50.0mm

04

939.640

64.0mm

05

939.780

78.0mm

06

These are multifaceted plates.
Can be used for Lisfranc,
Calcaneal cuboid, Naviculocunei-
form Fusion and more.

Screw used:
3.5mm LHS Screw,
3.5mm Cortex Screw,
2.7 Cortex Screws

MP Stick Plate



MP Mini Link Plate 2.4mm



Product Description

These are multifaceted plates
Can be used for Lisfranc,
Calcaneal cuboid, Naviculocunei-
form Fusion and more.

Screw used:
2.7mm LHS Screw,
2.4mm LHS Screw,
2.7mm Cortex Screw,
2.4mm Cortex Screw

MP Mini Link Plate 2.4mm

Codes	Range
973.002	02 Hole
973.003	03 Hole
973.004	04 Hole
973.006	06 Hole

MP Ankle Straight Plate 2.4mm



Product Description

MP Ankle Straight Plate is used for 1st Metatarsal Osteotomy for Metatarsal Deformity Correction, Tarsometatarsal Joint Arthrodesis, Medial Column Arthrodesis, Lateral Column Arthrodesis, Lisfranc Fracture/Dislocation, Metatarsal Fracture, Navicular Fracture, 5th Metatarsal Fracture.

Screw used: 2.7mm LHS Screw, 2.4mm LHS Screw

MP Ankle Straight Plate 2.4mm

Codes

978.006

MP Midfoot TMT Plate

Left

Right

Range

935.L.001

935.R.001

Small

935.L.002

935.R.002

Medium

Product Description

The Plate design and Screw placement allows very stiff fixation of Lisfranc fusion and other midfoot fixation

Screw used:
3.5mm LHS Screw
3.5mm Cortex Screw
2.7mm Cortex Screw

MP Midfoot TMT Plate



MP Mid-foot C Plate



Product Description

It is intended for fixation of osteotomies, fusions, non-unions, mal-unions, and fractures of the Mid-foot.

Screw used: 3.5mm LHS Screw, 3.5mm Cortex Screw, 5.0mm LHS Screw, 2.7mm Cortex Screw, 2.7 LHS Screw, 2.4mm Cortex Screw, 2.4mm LHS Screw.

MP Mid-foot C Plate

Left	Right	Range
934.L.002	934.R.002	Medium
934.L.003	934.R.003	Large

MP Universal X Plate

Product Description

It is used for treatment of Pilon Fractures, Trimalleolar Fractures, Bi-Malleolar Fracture, Cubiod Fracture, Navicular Fracture, Lisfranc, Syndesmosis Injury, Ankle Arthrodesis, Triple Arthrodesis, Calcaneal Cubiod, Metatarsocuneiform Fusion, Lapidus, Fibula Fractures, Talonavicularcuneiform Fusion, Radius Fracture, Ulna Fracture, Clavicle Fracture and more

Screw used: 3.5mm LHS Screw, 3.5mm Flat Head Cortex Screw

MP Universal X Plate



MP Mini Universal X Plate 2.4mm



Product Description

MP Mini Universal X Plate is intended for Multiplanar fixation.

Screw used:

2.7mm LHS Screw
2.4mm LHS Screw
2.7mm Cortex Screw
2.4mm Cortex Screw

MP Mini Universal X Plate 2.4mm

Codes

Range

974.001

Small

974.002

Medium

974.003

Large

MP First MTP Fusion Plate

Codes

Range

938.L.002

Left

938.R.002

Right

Product Description

It is used for 1st MTP Fusion and to treat hallux valgus that is caused by first TMT joint hypermobility.

Screw used: 3.5mm LHS Screw,
3.5mm Cortex Screw, 2.7mm
Cortex Screw.

MP First MTP Fusion Plate



MP Mini MTP Plate 2.4mm



Product Description

It is used for 1st MTP Fusion and to treat hallux valgus that is caused by first TMT joint hypermobility.

Screw used:

2.7mm LHS Screw, 2.4mm LHS Screw, 2.7mm Cortex Screw, 2.4mm Cortex Screw

MP Mini MTP Plate 2.4mm

Codes	Range
972.L.006	Left
972.R.006	Right

MP Mini Y Plate 2.4mm

Codes	Range
971.001	Small
971.002	Large

Product Description

It is intended to stabilize and aid in the repair of fractures, fusions, and osteotomies for small bones and bone fragments and also for Lapidus Arthrodesis,

Screw used: 2.7mm LHS Screw, 2.4mm LHS Screw, 2.7mm Cortex Screw, 2.4mm Cortex Screw

MP Mini Y Plate 2.4mm



MP Mini X Plate 2.4mm



Product Description

MP Mini X Plate Is intended for treatment of arthrodesis, proximal osteotomies and calcaneal osteotomies.

Screw used:
2.7mm LHS Screw
2.4mm LHS Screw

MP Mini X Plate 2.4mm

Codes	Range
975.004	04 Hole

MP Ankle Y Plate With Bar 3.5mm

Product Description

It provides fixation options for fractures of distal tibia and fibula. It can be used for minimally invasive Calcaneal Surgery

Screw used:
3.5mm LHS Screw
3.5mm Flat Head Cortex Screw

MP Ankle Y Plate With Bar 3.5mm



MP Ankle L Plate 2.4mm Obl Angle



Product Description

It addresses proximal and base metatarsal fractures and/or osteotomies, as well as Medial Column Arthrodesis,

Screw used:

2.7mm LHS Screw

2.4mm LHS Screw

MP Ankle L Plate 2.4mm Obl Angle

Codes	Range
977.L.007	Left
977.R.007	Right

MP Rhomboid Plate 3.5mm		
3.5mm Codes		
Left	Right	Range
979.L.001	979.R.001	Small
979.L.002	979.R.002	Medium
979.L.003	979.R.003	Large

Product Description

It is used for treatment of Lisfranc Arthrodesis, TN Arthrodesis, CC Arthrodesis and also for lateral column lengthening

Screw used:

3.5mm LHS Screw,

3.5mm Flat Head Cortex Screw

MP Rhomboid Plate 3.5mm



MP Ankle Arthrodesis Nail



Product Description

MP Ankle Arthrodesis Nail is indicated for tibiotalocalcaneal arthrodesis.

Screw used:
3.9 Locking Bolts
4.9mm Locking Bolts

MP Ankle Arthrodesis Nail

Codes			Range
9.0mm	10.00mm	11.00mm	
926.09.150	926.10.150	926.11.150	150
926.09.180	926.10.180	926.11.180	180
926.09.210	926.10.210	926.11.210	210
926.09.240	926.10.240	926.11.240	240
926.09.270	926.10.270	926.11.270	270
926.09.300	926.10.300	926.11.300	300
926.09.330	926.10.330	926.11.330	330

Cortex Flat Head Screw 3.5mm		Product Description
Codes	Length	Cortex Flat Head Screw 3.5mm is used to fix the Foot & Ankle plates.
925.010	10 mm	
925.012	12 mm	
925.014	14 mm	
925.016	16 mm	
925.018	18 mm	
925.020	20 mm	
925.022	22 mm	
925.024	24 mm	
925.026	26 mm	
925.028	28 mm	
925.030	30 mm	
925.032	32 mm	
925.034	34 mm	
925.036	36 mm	
925.038	38 mm	
925.040	40 mm	
925.042	42 mm	

Cortex Flat Head Screw 3.5mm



MP CC Screw 4.0mm



The MP CC Screw is used in common fractures, ankle arthrodesis, subtalar arthrodesis, and triple arthrodesis. Fractures of the fibula, medial malleolus,

Also used in treatment of tibial plateau and femoral intercondylar

Product Description

MP CC Screw 4.0mm Short Thread

4.0mm	Code	Range
TIT	610	10mm to 30mm (In steps of 02mm) 30mm to 70mm (In steps of 5mm)
SS	128	10mm to 30mm (In steps of 02mm) 30mm to 70mm (In steps of 5mm)

MP CC Screw 4.0mm Full Thread

4.0mm	Code	Range
TIT	611	10mm to 30mm (In steps of 02mm) 30mm to 70mm (In steps of 5mm)
SS	130	10mm to 30mm (In steps of 02mm) 30mm to 70mm (In steps of 5mm)

MP CC Screw 5.5mm			Product Description
Codes		Length	The MP CC Screw is used in common fractures, ankle arthrodesis, subtalar arthrodesis, double hindfoot arthrodesis and triple arthrodesis. Also used for fractures of the fibula, medial malleolus.. Also used in treatment of tibial plateau and femoral intercondylar
16mm	32mm		
940.025	-	25mm	
940.030	-	30mm	
940.035	941.035	35mm	
940.040	941.040	40mm	
940.045	941.045	45mm	
940.050	941.050	50mm	
940.055	941.055	55mm	
940.060	941.060	60mm	
940.065	941.065	65mm	
940.070	941.070	70mm	
940.075	941.075	75mm	
940.080	941.080	80mm	

MP CC Screw 5.5mm



MP CC Screw 7.0mm



Product Description

The MP Screw 7.0mm is commonly used in treatment of calcaneal fracture as well intended for fixation arthrodesis of the metatarsal-cuneiform, navicular-cuneiform, metatarsal-zuboid, talonavicular, and calcaneocuboid joints. It is also used for double hindfootarthrodesis and triple arthrodesis.

MP CC Screw 7.0mm

Codes		Length
16mm	32mm	
898.030	-	30mm
898.035	899.035	35mm
898.040	899.040	40mm
898.045	899.045	45mm
898.050	899.050	50mm
898.055	899.055	55mm
898.060	899.060	60mm
898.065	899.065	65mm
898.070	899.070	70mm
898.075	899.075	75mm
898.080	899.080	80mm
898.085	899.085	85mm
898.090	899.090	90mm
898.095	899.095	95mm
898.100	899.100	100mm
898.105	899.105	105mm
898.110	899.110	110mm
898.115	899.115	115mm
898.120	899.120	120mm
898.125	899.125	125mm
898.130	899.130	130mm



An innovative Plate System Designed for foot & ankle Surgery





Key Features

- Two Plane Locking at calcaneus (40% More stability)
- 6.8mm Compression at Tibiotalar Level
- Straight nail